

tactoset™

Bone Substitute Material

Caution: Federal law restricts this device to sale by or on the order of a physician

INSTRUCTIONS FOR USE

DESCRIPTION

Tactoset is a synthetic, biocompatible bone graft substitute material that hardens and converts to a poorly crystalline hydroxyapatite at body temperature. It is provided in single patient, single use kits in various volumes appropriate to the surgical site.

INTENDED USE / INDICATIONS

Tactoset is a synthetic, biocompatible bone graft substitute material that hardens and converts to a poorly crystalline hydroxyapatite at body temperature. It is indicated for filling bone voids or defects of the skeletal system (i.e. extremities and pelvis) that are not intrinsic to the stability of bony structure. These defects may be surgically created osseous defects or defects created from traumatic injury to the bone. The device provides an injectable, self-setting, osteoconductive bone graft substitute that resorbs and is replaced by the growth of new bone during the healing process.

DURATION OF ADMINISTRATION

Tactoset is intended for permanent implantation.

CONTRAINDICATIONS

Do not use this product if one or more of the following conditions are present:

- Existing acute or chronic infections, especially at the site of the operation.
- Nonviable bone.
- Areas where surrounding bone is not viable or not capable of supporting and anchoring the implant.
- Altered calcium metabolism.
- Metabolic bone disease.
- Immunologic abnormalities.
- Systemic disorders which result in poor wound healing.
- Inflammatory bone disease.
- Acute traumatic injuries with open wounds close to the defect which are likely to become infected.

METHOD OF STERILIZATION

- Powder is provided STERILE via gamma irradiation (Cobalt ⁶⁰Co).
- Liquid is provided STERILE via steam sterilization.
- Liquid syringe and packaging is provided STERILE via vaporized hydrogen peroxide.
- Contents are STERILE unless the barrier packaging is open or damaged; DO NOT USE if the package is open or damaged.
- Contents are non-pyrogenic.

PRODUCT STORAGE

Store at controlled room temperature within 15 to 25°C (59 to 77°F).

- The expiration date is printed on the outer package label and patient record labels.
- DO NOT USE expired product.

WARNINGS

- DISCARD UNUSED PORTIONS; Tactoset is a SINGLE PATIENT, SINGLE USE product.
- DO NOT RESTERILIZE; the safety and effectiveness of reused or resterilized Tactoset is unknown. Because Tactoset is indicated for use in defects that are not intrinsic to the stability of the bony structure, it is critical that adequate fixation be provided for unstable defects by other means.
- The safety and effectiveness for patients having received or to receive chemotherapy or radiation therapy at or near the implant site is not known.



tactoset™

injectable bone substitute

- The safety and effectiveness when used in conjunction with other legally marketed devices having similar indications is not known.
- The safety and effectiveness for use in children or elderly patients is not known.
- The effect in patients with documented renal disease is not known.
- The effect in patients with metabolic bone disease is not known.
- The effect in patients that are pregnant/nursing is not known.
- The effect in patients with cardiovascular disease precluding elective surgery is not known.
- The effect in patients having had infection during the last 3 months is not known.
- Care must be taken to prevent the creation of emboli. Highly pressurized application of Tactoset into a tightly confined space with ready venous or arterial access is not recommended, as the potential for formation of emboli is unknown.
- Prepare Tactoset using only the specified Mixing Solutions; the effect of preparing with other substances is unknown and may adversely affect product performance.

PRECAUTIONS

- Only for use by surgeons familiar with the material, appropriate surgical techniques, and bone repair procedures.
- Use aseptic technique to minimize the risk of infection.
- Mix with the specified volume of Mixing Solutions; deviations will alter the consistency of the material and may adversely affect the setting reaction and the effectiveness of the implant.
- Do not disturb the material after implantation as disruption may affect the characteristics of the hardened material.
- Do not irrigate the defect site immediately after implantation; wait until the material is hard to touch (about 10 minutes at 37 °C, 98 °F).
- Do not overfill the defect; Over-pressurizing the device may lead to extrusion beyond the site of intended application and damage to surrounding tissue. Remove any excess material within 2 minutes following implantation.
- Follow general surgical protocol regarding use of fixation.
- Post-operative use of a closed suction drain is recommended.
- Do not administer to patients with known hypersensitivity (allergy) to hyaluronate preparations.

POSSIBLE COMPLICATIONS

Re-operation to remove or replace the implant may be required occasionally due to medical reasons or device failure; if corrective action is not taken, one or more of the following complications may occur:

- Tissue thinning over implant site.
- Tenderness/redness/edema.
- Seroma/hematoma or infection.
- Swelling/fluid collection.
- Loss of contour.
- Migration, extrusion, dehiscence, fracture and sloughing of Tactoset can occur as a result of excessive trauma or post-operative load bearing.
- Neurovascular injuries due to surgical trauma.

MIXING INSTRUCTIONS

The kit components are identified as follows:

1. Mixing syringe containing Tactoset powder
2. Glass syringe containing Tactoset setting solution
3. Female-female Luer lock connector
4. 1 mL Luer lock syringes

STEP 1: HYDRATE TACTOSET POWDER

Hydration is achieved as follows:

- Transfer kit components into the sterile field using sterile technique.
- Discard packaging.
- Remove the cap from the glass syringe containing the setting solution.
- Remove the cap from the mixing syringe containing the powder (do not discard).
- Using a female-female Luer lock connector, attach the glass syringe containing the setting solution to the mixing syringe containing the powder.
- Hold the syringes vertically with the glass syringe containing the setting solution on top and inject the setting solution into the mixing syringe containing the powder.

- Remove the glass syringe and Luer lock connector and recap the mixing syringe.
- Remove the plunger rod sleeve (do not discard). Using the rod of the integrated mixing device, continuously mix the material from top to bottom of the mixing syringe for one minute.
- Reattach the plunger rod sleeve. Attach a Luer lock connector to the mixing syringe and expel residual air until a small drop of Tactoset is ejected.
- Attach a 1 mL Luer lock syringe. Fill the 1 mL syringe to 0.9 mL. Fill and repeat as needed with remaining syringes.

STEP 2: INJECT

















Working time is 7-18 minutes from when the powder and setting solution make contact; material must be implanted in this time.



Setting time is an inverse function of local body temperature; approximately 10 minutes at 37 °C (98 °F).

- Tactoset is intended to be injected. It should not be implanted as a putty.
- To inject, use an 11-15 Ga Delivery Cannula.
- Do not overfill the defect; remove excess material prior to hardening and avoid disturbing to allow proper setting.
- Allow approximately 10 minutes after injection for Tactoset to harden prior to irrigating.
- Close the defect and follow general surgical protocol regarding post-operative use of closed suction drain.

	Consult Instructions for Use		Do Not Use If Package Is Damaged
	Caution		Lot Number
	Do Not Reuse		Expiration Date
	Sterilized using Radiation		Manufacturer
	Sterilized using Vaporized Hydrogen Peroxide		Keep Away From Sunlight
	Sterilized using Steam		Keep Dry
	CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician		Temperature Limitation

Anika Therapeutics, Inc.
 32 Wiggins Avenue
 Bedford, MA 01730
 USA
 Phone: (888) 721-1600
 customerservice@anikatherapeutics.com
<http://www.anikatherapeutics.com/>