GFS Depth Gauge
Medidor de profundidad GFS
Jauge de profondeur GFS
GFS-Tiefenlehre
Misuratore di profondità GFS
Medidor de Profundidade GFS
Μετρητής Βάθους GFS
GFS gylio matuoklis
GFS dybdemåler
GFS djupmätare
Jožă de adâncime GFS
Głębokościomierz GFS
Манометър за дълбочина GFS

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Wichtige Produktinformationen
Informazioni importanti sul prodotto
Informações Importantes sobre o Produto
Σημαντικές πληροφορίες για το προϊόν
Svarbi informacija apie gaminį
Viktig produktinformasjon
Viktig produktinformation
Informații referitoare la produs
Ważne informacje o produkcie
Важна информация за продукта

Directions for Use
Instrucciones de uso
Instructions d'utilisation
Gebrauchsanweisung
Istruzioni per l'uso
Instruções para o uso
Οδηγίες χρήσης
Naudojimo instrukcija
Bruksanvisning
Anvisningar för användning
Instrucțiuni de utilizare
Wskazówki dotyczące użytkowania
Инструкции за употреба

30389 Rev 6
Eff. 07/24

 Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

 **ANIKA**  2797
www.anika.com
(941)755-7965

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GFS Depth Gauge (EN)

1. Intended Use:

The GFS Depth Gauge is indicated for use with the Anika Graft Fixation System (GFS) and may be used to measure the overall length of a tunnel as well as for measuring and marking the graft.

2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. Do not use this product if it is bent or there is visible damage to the shaft or laser markings.

3. Material Specifications:

The GFS Depth Gauge is manufactured from a Stainless-Steel Alloy. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

5. Cleaning:

- A. Completely submerge used instruments in enzyme solution and allow to soak for 20 minutes. Use a soft bristled, nylon brush to gently scrub the instrument until all visible soil has been removed. Particular attention must be given to crevices, lumen, mated surfaces, connections and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft bristled brush (i.e. pipe cleaner brush).
- B. Remove instrument from enzymatic solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas .
- C. Place prepared cleaning agents in a sonication unit. Completely submerge instrument in cleaning solution and sonicate for 10 minutes at 45-50 kHz.
- D. Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the instrument or in the rinse stream. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- E. Repeat steps A through D a second time and ensure the instrument is free of blood or soil before moving to step F.
- F. Remove excess moisture from the instrument with a clean, absorbent, and non-shedding wipe.

6. Sterilization:

This device is not provided sterile. Devices shall be double wrapped in wrapping material or pouches, double wrapped while in a container or tray, or placed within a reusable sterilization container. Combining two or more systems may hinder the efficacy of the sterilization to an unacceptable level. Reusable devices shall be sterilized with a load size which does not exceed 10,630 grams. The following table provides recommended minimum sterilization parameters that have been validated by Anika Therapeutics to provide a 10⁻⁶ sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes ¹
	134°C / 273°F	3 minutes	
	132°C / 270°F	18 minutes ²	50 minutes ¹


¹ Drying times vary according to load size and should be increased for larger loads.

² This cycle is not for use in the United States. These are disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern regarding TSE/CDJ contamination. This cycle is not to be used for the inactivation of TSE/CJD contamination.

7. Directions for Use (Soft Tissue Graft):

- A. To measure the overall length of a tunnel, place GFS Depth Gauge in the tunnel and insert the hooked end of the gauge until it catches on the cortical bone at the distal end of the tunnel. Read the graduations on the GFS Depth Gauge to determine the length of the tunnel.
- B. To aid in marking the graft, assemble the GFS/Graft construct per GFS Directions for Use and orient it in such a way that the axis of the titanium implant is parallel to the graft. Place the wedge shaped end of the GFS at the zero point on the GFS Depth Gauge and make a mark on the graft with a surgical marker that corresponds to the overall length of the femoral

tunnel. This will indicate the point at which the wedge shaped end of the GFS has passed the lateral femoral cortex allowing the implant to flip into position. See Drawing below.

8. Adverse/ Reportable Event:

Any serious incident that occurs in relation to the device should be reported to both Anika Therapeutics and the applicable authority using the following contact information.

Country	Authority	Email	Website
All	Anika Therapeutics	ParcusCS@Anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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