Anika Tendon Harvester
Extractor de tendones Anika
Récupérateur de tendon Anika
Anika Sehenschneider
Raccoglitore di tendini Anika
Colhedora de Tendão Anika
Συλλέκτης Τενόντων Anika
„Anika“ sausgyslių ištraukimo prietaisas
Anika senestripper
Anika insamlingsverktug av sena
Colector de tendon Anika
Kombajn do ściągacien Anika
Уред за извличане на сухожилия Anika

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Wichtige Produktinformationen
Informazioni importanti sul prodotto
Informações Importantes sobre o Produto
Σημαντικές πληροφορίες για το προϊόν
Svarbi informacija apie gaminį
Viktig produktinformasjon
Viktig produktinformation
Informații referitoare la produs
Ważne informacje o produkcie
Важна информация за продукта

Directions for Use
Instrucciones de uso
Instructions d'utilisation
Gebrauchsanweisung
Istruzioni per l'uso
Instruções para o uso
Οδηγίες χρήσης
Naudojimo instrukcija
Bruksanvisning
Anvisningar för användning
Instrucțiuni de utilizare
Wskazówki dotyczące użytkowania
Инструкции за употреба

30851 Rev 6
Eff. 08/24

 Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

 ANIKA www.anika.com
(941)755-7965

 2797
EC REP
MDSS GmbH
Schiffgraben 41
30175 Hanover, Germany

Anika Tendon Harvester (EN)

1. Intended Use:

The Anika Tendon Harvester is used to remove tendons (e.g. semitendinosus and gracilis) for use in repairing or replacing damaged tendons and ligaments.

2. Warnings:

- a. This product is intended for use by or on the order of a physician.
- b. No component of this instrument is intended to be left in the body.
- c. Never use this device if it has been bent or damaged in any way.
- d. Caution: Common risks consistent with surgery may include: bone/tissue damage, revision surgery, soft tissue irritation, infection, allergic reactions.

3. Material:

This device is manufactured from stainless steel and aluminum alloys. The materials used in the manufacture of this device are radiopaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

Instruments intended to be reused should be inspected prior to each use to assess damage and wear. Signs of damage and wear include surface discoloration, surface corrosion, surface damage, and general bending, cracking, and fractures. Cutting, drilling, and other instruments intended to have sharp edges will become dull over time, diminishing their ability to function as designed. Review sharp instruments for dullness or damage at the cutting/sharp feature. If the instrument does not pass inspection, do not use the instrument, and notify an Anika representative.

5. Storage:

Products must be stored in a dry place.

6. Cleaning:

- a. Completely submerge used instruments in enzyme solution and allow to soak for 20 minutes. Use a soft bristled, nylon brush to gently scrub the instrument until all visible soil has been removed. Particular attention must be given to crevices, lumen, mated surfaces, connections and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft bristled brush (i.e. pipe cleaner brush).
- b. Remove instrument from enzymatic solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- c. Place prepared cleaning agents in a sonication unit. Completely submerge instrument in cleaning solution and sonicate for 10 minutes at 45-50 kHz.
- d. Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the instrument or in the rinse stream. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- e. Repeat steps A through D a second time and ensure the instrument is free of blood or soil before moving to step F.
- f. Remove excess moisture from the instrument with a clean, absorbent, and non-shedding wipe.

7. Sterilization:

This device is not provided sterile. Devices shall be double wrapped in wrapping material or pouches, double wrapped while in a container or tray, or placed within a reusable sterilization container. Combining two or more systems may hinder the efficacy of the sterilization to an unacceptable level. Reusable devices shall be sterilized with a load size which does not exceed 10,630 grams. The following table provides recommended minimum sterilization parameters that have been validated by Anika Therapeutics to provide a 10^{-6} sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes ¹
	134°C / 273°F	3 minutes	

	132°C / 270°F	18 minutes ²	50 minutes ¹
--	---------------	-------------------------	-------------------------

¹ Drying times vary according to load size and should be increased for larger loads.

² This cycle is not for use in the United States. These are disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern regarding TSE/CDJ contamination. This cycle is not to be used for the inactivation of TSE/CJD contamination.

After sterilization aseptically transfer device into the sterile field in accordance with the institution's policy.

8. Directions for use:

- Aseptically transfer device into sterile field and identify the tendon to be harvested taking care to select a harvester best matched to the tendon diameter.
- Anika Tendon Harvesters come in two versions, closed and slotted (open). When using the closed harvester, it is advisable to place sutures through the tendon, at its insertion. Sutures will aid in providing the counter-traction necessary to harvest the tendon.
- The accessory insertion of the semitendinosus tendon should be identified and transected to avoid tendon damage at harvest.

Note: Knee flexion may reduce the risk of injury to the saphenous nerve as it crosses the gracilis tendon.

- Once the tendon is released, remove the instrument and tendon from the wound and prepare the graft according to the prescribed surgical technique.

9. Device Disposal:

Devices that have evidence of damage must be properly disposed of in accordance with the institution's policy.

10. Adverse/ Reportable Event:

Any serious incident that occurs in relation to the device should be reported to both Anika Therapeutics and the applicable authority using the following contact information.

Country	Authority	Email	Website
All	Anika Therapeutics	ParcusCS@Anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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